

KERN 346-83



Individual weight, OIML Class M1, 500 kg, Heavy Duty Weight, cast iron lacquered



Category

| | |
|------------------|-------------------|
| Brand | KERN |
| Product category | Individual weight |
| Product group | OIML M1 |
| Product family | 346 |

Measuring System

| | |
|---|--------|
| Test weight - nominal value | 500 kg |
| Test weight - volume (cm ³) | 69444 |
| Test weight - volume uncertainty (cm ³) | 5787 |
| Permissible tolerance according to OIML (+/-) | 25 g |
| Test weight - Density [Min] (OIML) | 3000 |
| Weight OIML class | M1 |

Construction

| | |
|---------------------------|---------------------|
| Shape weight | Heavy Duty Weight |
| Dimensions weight (W×D×H) | 750×500×314 mm |
| Material | cast iron lacquered |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 2 d |
| Dimensions packaging (W×D×H) | 750×500×314 mm |
| Net weight | 500 kg |
| Shipping method | Freight |
| Net weight approx. | 500 kg |
| Gross weight approx. | 500 kg |
| Shipping weight | 500 kg |

Services (optional)

| | |
|--------------------------------------|---------|
| Article number for DAkkS calibration | 962-693 |
|--------------------------------------|---------|

Pictograms

STANDARD

