

KERN 346-84



Individual weight, OIML Class M1, 1000 kg, Heavy Duty Weight, cast iron lacquered



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

Measuring System

Test weight - nominal value	1000 kg
Test weight - volume (cm ³)	138889
Test weight - volume uncertainty (cm ³)	11574
Permissible tolerance according to OIML (+/-)	50 g
Test weight - Density [Min] (OIML)	3000
Weight OIML class	M1

Construction

Shape weight	Heavy Duty Weight
Dimensions weight (W×D×H)	750×500×500 mm
Material	cast iron lacquered

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	750×500×500 mm
Net weight	1000 kg
Shipping method	Freight
Net weight approx.	1000 kg
Gross weight approx.	1000 kg
Shipping weight	1000 kg

Services (optional)

Article number for DAkkS calibration	962-694
--------------------------------------	---------

Pictograms

STANDARD

