

KERN 346-85



Individual weight, OIML Class M1, 2000 kg, Heavy Duty Weight, cast iron lacquered



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

Measuring System

Test weight - nominal value	2000 kg
Test weight - volume (cm ³)	277778
Test weight - volume uncertainty (cm ³)	23148
Permissible tolerance according to OIML (+/-)	100 g
Test weight - Density [Min] (OIML)	3000
Weight OIML class	M1

Construction

Shape weight	Heavy Duty Weight
Dimensions weight (W×D×H)	1000×750×500 mm
Material	cast iron lacquered

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	1000×750×500 mm
Net weight	2000 kg
Shipping method	Freight
Net weight approx.	2000 kg
Gross weight approx.	2000 kg
Shipping weight	2000 kg

Services (optional)

Article number for DAkkS calibration	962-695
--------------------------------------	---------

Pictograms

STANDARD

