

# KERN 347-009-400

**KERN**

Plastic box for individual weights (Cylindrical/compact, 1 mg - 500 mg (einzel))



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Plastic box

## Measuring System

Case suitable for OIML classes	E1 - M2
Case suitable for nominal weight values	1 mg - 500 mg (einzel)
Suitable weight shape for case	Cylindrical/compact

## Construction

Material	plastic
Dimension (Ø×H)	24×51 mm

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	24×24×51 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Pictograms

STANDARD

