


# KERN 347-022-005

# KERN

Weight container for test weights, calibrated to OIML class M1



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Weight container
Product family	347

## Measuring System

Test weight - nominal value	20 kg
Permissible tolerance according to OIML (+/-)	1 mg
Weight OIML class	M1

## Construction

Shape weight	weight container
Max. total weight	160 kg
Material	stainless steel

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	1000×500×500 mm
Net weight	20 kg
Shipping method	Parcel service
Net weight approx.	20 kg
Gross weight approx.	20 kg
Shipping weight	50 kg

## Services (optional)

Article number for DAkkS calibration	962-645
--------------------------------------	---------