


KERN 347-052-010

KERN

Weight container for test weights, calibrated to OIML class M1



Pictograms

STANDARD



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Weight container
Product family	347

Measuring System

Test weight - nominal value	50 kg
Permissible tolerance according to OIML (+/-)	2,5 mg
Weight OIML class	M1

Construction

Shape weight	weight container
Material	stainless steel

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	1000×500×500 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg

Services (optional)

Article number for DAkkS calibration	962-646
--------------------------------------	---------