

# KERN 347-06

**KERN**

Individual weight, OIML Class M1, 50 g, Knob, stainless steelfine turned (OIML)



## Services (optional)

Article number for conformity assessment	952-636
Article number for DAkkS calibration	962-636

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

## Measuring System

Test weight - nominal value	50 g
Test weight - volume (cm <sup>3</sup> )	6,25
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0781
Permissible tolerance according to OIML (+/-)	3 mg
Test weight - Density [Min] (OIML)	4000
Weight OIML class	M1

## Construction

Shape weight	Knob
Dimensions (Ø×H)	18×29,5 mm
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	18×18×29 mm
Net weight	0,05 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg