

KERN 347-07

KERN

Individual weight, OIML Class M1, 100 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for conformity assessment	952-637
Article number for DAkkS calibration	962-637

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - nominal value	100 g
Test weight - volume (cm ³)	12,5
Test weight - volume uncertainty (cm ³)	0,1563
Permissible tolerance according to OIML (+/-)	5 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

Construction

Shape weight	Knob
Dimensions (Ø×H)	22×39 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	22×22×39 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg