

KERN 347-08

Individual weight, OIML Class M1, 200 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for conformity assessment	952-638
Article number for DAkkS calibration	962-638

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - nominal value	200 g
Test weight - volume (cm ³)	25
Test weight - volume uncertainty (cm ³)	0,3125
Permissible tolerance according to OIML (+/-)	10 mg
Test weight - Density [Min] (OIML)	4000
Weight OIML class	M1

Construction

Shape weight	Knob
Dimensions (Ø×H)	28×48,5 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	28×28×48 mm
Net weight	0,2 kg
Shipping method	Parcel service
Net weight approx.	0,20 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg