

# KERN 347-09

Individual weight, OIML Class M1, 500 g, Knob, stainless steelfine turned (OIML)



## Services (optional)

Article number for conformity assessment	952-639
Article number for DAkkS calibration	962-639

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

## Measuring System

Test weight - nominal value	500 g
Test weight - volume (cm <sup>3</sup> )	62,5
Test weight - volume uncertainty (cm <sup>3</sup> )	0,7813
Permissible tolerance according to OIML (+/-)	25 mg
Test weight - Density [Min] (OIML)	4000
Weight OIML class	M1

## Construction

Shape weight	Knob
Dimensions (Ø×H)	38×66 mm
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	38×38×66 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg