

KERN 347-130-400

KERN

Plastic box for individual weights (Cylindrical/compact, 5 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Plastic box

Measuring System

Case suitable for OIML classes	F1 - M3
Case suitable for nominal weight values	5 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Material	plastic
Dimension (Ø×H)	100×165 mm

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	100×100×165 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,3 kg

Pictograms

STANDARD

