

KERN 347-141

KERN

Individual weight, OIML Class M1, 10 kg, Check weight, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-644

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - nominal value	10 kg
Test weight - volume (cm ³)	1250
Test weight - volume uncertainty (cm ³)	15,625
Permissible tolerance according to OIML (+/-)	500 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	M1

Construction

Shape weight	Check weight
Dimensions (Ø×H)	136×132 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	137×137×132 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg