

KERN 347-151

KERN

Individual weight, OIML Class M1, 20 kg, Check weight, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-645

Pictograms

STANDARD



Category

| | |
|------------------|-------------------|
| Brand | KERN |
| Product category | Individual weight |
| Product group | OIML M1 |
| Product family | 347 |

Measuring System

| | |
|---|-------|
| Test weight - nominal value | 20 kg |
| Test weight - volume (cm ³) | 2500 |
| Test weight - volume uncertainty (cm ³) | 31,25 |
| Permissible tolerance according to OIML (+/-) | 1 g |
| Test weight - Density [Min] (OIML) | 7390 |
| Test weight - Density [Max] (OIML) | 8730 |
| Weight OIML class | M1 |

Construction

| | |
|------------------|-----------------------------------|
| Shape weight | Check weight |
| Dimensions (Ø×H) | 136×217 mm |
| Material | stainless steelfine turned (OIML) |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 137×137×217 mm |
| Net weight | 20 kg |
| Shipping method | Parcel service |
| Net weight approx. | 20 kg |
| Gross weight approx. | 20 kg |
| Shipping weight | 20 kg |