

KERN 347-161

KERN

Individual weight, OIML Class M1, 50 kg, Check weight, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-646

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - nominal value	50 kg
Test weight - volume (cm ³)	6250
Test weight - volume uncertainty (cm ³)	78,125
Permissible tolerance according to OIML (+/-)	2,5 g
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	M1

Construction

Shape weight	Check weight
Dimensions (Ø×H)	198×249,5 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	198×198×249 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg