

KERN 347-445-100

Supporting bar, OIML Class M1, 10 g, aluminium

KERN



Services (optional)

Article number for DAkkS calibration	962-634
--------------------------------------	---------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	Supporting Bar
Product family	347

Measuring System

Test weight - nominal value	10 g
Test weight - volume (cm ³)	1,1905
Test weight - volume uncertainty (cm ³)	0,0241
Permissible tolerance according to OIML (+/-)	2 mg
Weight OIML class	M1

Construction

Shape weight	Supporting bar (Type A)
Max. total weight	0,2 kg
Max. single weight	0,1 kg
Dimensions (Ø×H)	13×113,5 mm
Material	aluminium

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	13×13×114 mm
Net weight	0,01 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg