

KERN 347-5 15-400

Plastic box for individual weights (Cylindrical/compact, 500 g - 1 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Plastic box

Measuring System

Case suitable for OIML classes	M1 - M3
Case suitable for nominal weight values	500 g - 1 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Material	plastic
Dimension (Ø×H)	80×600 mm

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	80×80×600 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,8 kg

Pictograms

STANDARD

