

KERN 348-226

KERN

Set of milligram weights, M1, 1 mg - 500 mg, Platelet, stainless steel polished (OIML), in aluminium case

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	348

Measuring System

Test weight - nominal value	1 mg - 500 mg
Weight OIML class	M1

Construction

Shape weight	Platelet
Dimension (W×D×H)	165×158×37 mm
Case of weight set	in aluminium case
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	165×158×37 mm
Net weight	0,32 kg
Shipping method	Parcel service
Net weight approx.	0,35 kg
Gross weight approx.	0,40 kg
Shipping weight	0,39 kg

Services (optional)

Article number for conformity assessment	952-650
Article number for DAkkS calibration	962-650