


KERN 353-420-210

Cover for box of microgramme weights

KERN



Category

Brand	KERN
Product category	Test weights accessories
Product group	Case cover
Product family	353

Measuring System

Case suitable for nominal weight values	1 mg - 500 mg
Suitable weight shape for case	Cylindrical/compact

Construction

Material	coated metal
----------	--------------

Packing & Shipping

Readability force [d] (N)	1 d
Net weight	0,05 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,05 kg

Pictograms