

KERN 353-460-200

Wooden box for weight sets (Cylindrical/compact, 1 mg - 1 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F2 + M1 + M2
Case suitable for nominal weight values	1 mg - 1 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	230×180×112 mm
Material	beech

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	230×180×112 mm
Net weight	2 kg
Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,6 kg
Shipping weight	2,5 kg

Pictograms

STANDARD

