

# KERN 354-02

**KERN**

Weight set, M2, 1 g - 50 g, Knob, stainless steel fine turned, in wooden case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML M2
Product family	354

## Measuring System

Test weight - nominal value	1 g - 50 g
Weight OIML class	M2

## Construction

Shape weight	Knob
Dimension (WxDxH)	129x86x40 mm
Number of weights	8
Case of weight set	in wooden case
Material	stainless steel fine turned

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	129x86x40 mm
Net weight	0,26 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,30 kg
Shipping weight	0,26 kg

## Services (optional)

Article number for DAkkS calibration	962-615
--------------------------------------	---------