

# KERN 354-03

**KERN**

Weight set, M2, 1 g - 100 g, Knob, stainless steelfine turned (OIML), in wooden case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML M2
Product family	354

## Measuring System

Test weight - nominal value	1 g - 100 g
Weight OIML class	M2

## Construction

Shape weight	Knob
Dimension (W×D×H)	151×91×49 mm
Number of weights	9
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	151×91×49 mm
Net weight	0,51 kg
Shipping method	Parcel service
Net weight approx.	0,55 kg
Gross weight approx.	0,55 kg
Shipping weight	0,51 kg

## Services (optional)

Article number for DAkkS calibration	962-616
--------------------------------------	---------