

KERN 354-036

KERN

Weight set, M2, 1 g - 100 g, Knob, stainless steel fine turned, in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M2
Product family	354

Measuring System

Test weight - nominal value	1 g - 100 g
Weight OIML class	M2

Construction

Shape weight	Knob
Dimension (W×D×H)	165×125×75 mm
Case of weight set	in aluminium case
Material	stainless steel fine turned

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	165×125×75 mm
Net weight	0,65 kg
Shipping method	Parcel service
Net weight approx.	0,65 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg

Services (optional)

Article number for DAkkS calibration	962-616
--------------------------------------	---------