

KERN 354-05

KERN

Weight set, M2, 1 g - 500 g, Knob, stainless steelfine turned (OIML), in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M2
Product family	354

Measuring System

Test weight - nominal value	1 g - 500 g
Weight OIML class	M2

Construction

Shape weight	Knob
Dimension (W×D×H)	168×98×78 mm
Number of weights	12
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	168×98×78 mm
Net weight	1,71 kg
Shipping method	Parcel service
Net weight approx.	1,8 kg
Gross weight approx.	1,8 kg
Shipping weight	1,71 kg

Services (optional)

Article number for DAkkS calibration	962-618
--------------------------------------	---------