

KERN 354-07

KERN

Weight set, M2, 1 g - 2 kg, Knob, stainless steelfine turned (OIML), in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M2
Product family	354

Measuring System

Test weight - nominal value	1 g - 2 kg
Weight OIML class	M2

Construction

Shape weight	Knob
Dimension (WxDxH)	305x180x140 mm
Number of weights	15
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	305x180x140 mm
Net weight	9,01 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	9,01 kg

Services (optional)

Article number for DAkkS calibration	962-620
--------------------------------------	---------