

# KERN 354-430-200

Wooden box for weight sets (Cylindrical/compact, 1 g - 100 g)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	F2 + M1 + M2 + M3
Case suitable for nominal weight values	1 g - 100 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	154×61×47 mm
Material	beech

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	155×60×45 mm
Net weight	0,3 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,30 kg
Shipping weight	0,3 kg

## Pictograms

STANDARD

