

KERN 354-440-200

Wooden box for weight sets (Cylindrical/compact, 1 g - 200 g)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F2 + M1 + M2 + M3
Case suitable for nominal weight values	1 g - 200 g
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	140×80×70 mm
Material	beech

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	125×80×70 mm
Net weight	0,43 kg
Shipping method	Parcel service
Net weight approx.	0,45 kg
Gross weight approx.	0,45 kg
Shipping weight	0,43 kg

Pictograms

STANDARD

