

# KERN 362-96

**KERN**

Weight set, M3, 1 g - 1 kg, stainless steel (1 g - 50 g) and cast iron lacquered ( $\geq 100$  g), button ( $\leq 50$  g,  $\geq 500$  g) and cylinders (100 g, 200 g), in wood block



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML M3
Product family	362

## Measuring System

Test weight - nominal value	1 g - 1 kg
Weight OIML class	M3

## Construction

Shape weight	Knob
Dimension (WxDxH)	260x95x50 mm
Number of weights	13
Case of weight set	in wood block
Material	cast iron lacquered

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	260x95x50 mm
Net weight	2,75 kg
Shipping method	Parcel service
Net weight approx.	2,8 kg
Gross weight approx.	2,8 kg
Shipping weight	2,75 kg

## Services (optional)

Article number for DAkkS calibration	962-619
--------------------------------------	---------