

KERN 362-970-300

KERN

Wooden box for weight sets (Cylindrical/compact, 1 g - 2 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	M3
Case suitable for nominal weight values	1 g - 2 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	340×120×60 mm
Material	beech

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	340×120×60 mm
Net weight	1,3 kg
Shipping method	Parcel service
Net weight approx.	1,4 kg
Gross weight approx.	1,4 kg
Shipping weight	1,3 kg

Pictograms

STANDARD

