

KERN 362-99

KERN

Weight set, M3, 1 g - 10 kg, stainless steel (1 g - 50 g) and cast iron lacquered (≥ 100 g), button (≤ 50 g, ≥ 500 g) and cylinders (100 g, 200 g), in wood block



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M3
Product family	362

Measuring System

Test weight - nominal value	1 g - 10 kg
Weight OIML class	M3

Construction

Shape weight	Knob
Dimension (W×D×H)	510×150×70 mm
Number of weights	17
Case of weight set	in wood block
Material	cast iron lacquered

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	510×150×70 mm
Net weight	23,5 kg
Shipping method	Parcel service
Net weight approx.	24 kg
Gross weight approx.	24 kg
Shipping weight	23,5 kg

Services (optional)

Article number for DAkkS calibration	962-622
--------------------------------------	---------