

KERN 366-91

KERN

Individual weight, OIML Class M3, 100 g, Cylinders, cast iron lacquered



Services (optional)

Article number for DAkkS calibration 962-637

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	366

Measuring System

Test weight - nominal value	100 g
Test weight - volume (cm ³)	14,0845
Test weight - volume uncertainty (cm ³)	1,1902
Permissible tolerance according to OIML (+/-)	50 mg
Test weight - Density [Min] (OIML)	1500
Weight OIML class	M3

Construction

Shape weight	Cylinders
Dimensions (Ø×H)	31×22 mm
Material	cast iron lacquered

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	31×31×22 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg