

# KERN 366-92

**KERN**

Individual weight, OIML Class M3, 200 g, Cylinders, cast iron lacquered



## Services (optional)

Article number for DAkkS calibration 962-638

## Pictograms

### STANDARD



## Category

|                  |                   |
|------------------|-------------------|
| Brand            | KERN              |
| Product category | Individual weight |
| Product group    | OIML M3           |
| Product family   | 366               |

## Measuring System

|   |        |
|---|--------|
| Test weight - nominal value                         | 200 g  |
| Test weight - volume (cm <sup>3</sup> )             | 28,169 |
| Test weight - volume uncertainty (cm <sup>3</sup> ) | 2,3805 |
| Permissible tolerance according to OIML (+/-)       | 100 mg |
| Test weight - Density [Min] (OIML)                  | 1500   |
| Weight OIML class                                   | M3     |

## Construction

|                  |                     |
|------------------|---------------------|
| Shape weight     | Cylinders           |
| Dimensions (Ø×H) | 39×28 mm            |
| Material         | cast iron lacquered |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Readability force [d] (N)    | 1 d            |
| Dimensions packaging (W×D×H) | 39×39×28 mm    |
| Net weight                   | 0,2 kg         |
| Shipping method              | Parcel service |
| Net weight approx.           | 0,20 kg        |
| Gross weight approx.         | 0,20 kg        |
| Shipping weight              | 0,2 kg         |