

# KERN 366-93

**KERN**

Individual weight, OIML Class M3, 500 g, Knob, cast iron lacquered



## Services (optional)

Article number for DAkkS calibration 962-639

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	366

## Measuring System

Test weight - nominal value	500 g
Test weight - volume (cm <sup>3</sup> )	70,4225
Test weight - volume uncertainty (cm <sup>3</sup> )	5,9512
Permissible tolerance according to OIML (+/-)	250 mg
Test weight - Density [Min] (OIML)	1500
Weight OIML class	M3

## Construction

Shape weight	Knob
Dimensions (Ø×H)	46×68 mm
Material	cast iron lacquered

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	46×46×68 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,50 kg
Shipping weight	0,5 kg