

KERN 366-95

KERN

Individual weight, OIML Class M3, 2 kg, Knob, cast iron lacquered



Services (optional)

Article number for DAkkS calibration 962-642

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	366

Measuring System

Test weight - nominal value	2 kg
Test weight - volume (cm ³)	281,6901
Test weight - volume uncertainty (cm ³)	23,8048
Permissible tolerance according to OIML (+/-)	1 g
Test weight - Density [Min] (OIML)	1500
Weight OIML class	M3

Construction

Shape weight	Knob
Dimensions (Ø×H)	70×113 mm
Material	cast iron lacquered

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	70×70×113 mm
Net weight	2 kg
Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg