

KERN 366-96

KERN

Individual weight, OIML Class M3, 5 kg, Knob, cast iron lacquered



Services (optional)

Article number for DAkkS calibration 962-643

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	366

Measuring System

Test weight - nominal value	5 kg
Test weight - volume (cm ³)	704,2254
Test weight - volume uncertainty (cm ³)	59,512
Permissible tolerance according to OIML (+/-)	2,5 g
Test weight - Density [Min] (OIML)	1500
Weight OIML class	M3

Construction

Shape weight	Knob
Dimensions (Ø×H)	94×146 mm
Material	cast iron lacquered

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	94×94×146 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg