

# KERN 366-97

**KERN**

Individual weight, OIML Class M3, 10 kg, Knob, cast iron lacquered



## Services (optional)

Article number for DAkkS calibration 962-644

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	366

## Measuring System

Test weight - nominal value	10 kg
Test weight - volume (cm <sup>3</sup> )	1408,4507
Test weight - volume uncertainty (cm <sup>3</sup> )	119,024
Permissible tolerance according to OIML (+/-)	5 g
Test weight - Density [Min] (OIML)	1500
Weight OIML class	M3

## Construction

Shape weight	Knob
Dimensions (Ø×H)	121×193 mm
Material	cast iron lacquered

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	121×121×193 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg