

KERN ALT-A02

KERN

Set for density determination of liquids and solids with density $\leq/\geq 1$, only for models with readout [d] = 0,001 g



Category

Brand	KERN
Product category	Balance

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	390×290×135 mm
Net weight	0,66 kg
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	1,6 kg
Shipping weight	3,1 kg

Pictograms

STANDARD

