

# KERN ATB-US02

Ultrasound probe; 5 MHz,  $\varnothing$  12 mm

**KERN**



## Category

Brand	Sauter
Product category	Measurementtechnology
Product group	Probe

## Construction

Material	aluminium
----------	-----------

## Environmental conditions

Storage temperature [Min]	-10 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	150×100×2 mm
Net weight	0,101 kg
Shipping method	Parcel service
Net weight approx.	0,15 kg
Gross weight approx.	0,15 kg
Shipping weight	0,104 kg

## Pictograms

### STANDARD

