

# KERN ATB-US02

Ultrasone sonde; 5 MHz,  $\varnothing$  12 mm

**KERN**



## Categorie

Merk	Sauter
Productcategorie	Meettechniek
Productgroep	Sonde

## Ontwerp

Materiaal	aluminium
-----------	-----------

## Milieuomstandigheden

Opslagtemperatuur [Min]	-10 °C
Opslagtemperatuur [Max]	40 °C

## Verpakking & verzending

Afreesbaarheid kracht [d] (N)	1 d
Afmetingen verpakking (B×D×H)	150×100×2 mm
Nettogewicht	0,101 kg
Verzendmethode	Pakketdienst
Nettogewicht ca.	0,15 kg
Brutogewicht ca.	0,15 kg
Verzendgewicht	0,104 kg

## Pictogrammen

### STANDARD

