

KERN ATB-US06

Ultrasound probe; 5 MHz, \varnothing 8 mm

KERN



Category

Brand	Sauter
Product category	Measurementtechnology
Product group	Probe

Construction

Material	aluminium
----------	-----------

Environmental conditions

Storage temperature [Min]	-10 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	100×80×2 mm
Net weight	0,044 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Pictograms

STANDARD

