

KERN ATU-US01

Ultrasound probe; 2,5 MHz, \varnothing 14 mm

KERN



Category

Brand	Sauter
Product category	Measurementtechnology
Product group	Probe

Construction

Dimensions completely mounted (W×D×H)	57×35×35 mm
Material	aluminium

Environmental conditions

Storage temperature [Min]	-10 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	150×100×20 mm
Net weight	0,0555 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

Pictograms

STANDARD

