

# KERN ATU-US09

**KERN**

Ultrasound probe; 5 MHz,  $\varnothing$  10 mm



## Category

Brand	Sauter
Product category	Measurementtechnology
Product group	Probe

## Construction

Dimensions completely mounted (W×D×H)	57×35×35 mm
Material	aluminium

## Environmental conditions

Storage temperature [Min]	-10 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	150×210×10 mm
Net weight	0,05 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Pictograms

STANDARD

