

KERN ATU-US12

Ultrasound probe; 5 MHz, \varnothing 10 mm

KERN



Category

Brand	Sauter
Product category	Measurementtechnology
Product group	Probe

Construction

Design	round
Dimensions completely mounted (WxDxH)	10x10x25 mm
Material	aluminium

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	70x100x2 mm
Net weight	0,068 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

Pictograms

STANDARD

