

KERN MLB-A01A

KERN

Sample plates aluminium, \varnothing 90 mm, unit of 80 pieces



Category

Brand	KERN
Product category	Balance
Product group	Bowl

Environmental conditions

Storage temperature [Min]	-100 °C
Storage temperature [Max]	100 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	105×105×80 mm
Net weight	0,3 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,35 kg
Shipping weight	0,34 kg

Pictograms

STANDARD

