

KERN MLB-A 11

KERN

Temperature calibration set consists of measuring sensor and display device



Category

Brand	KERN
Product category	Balance
Product group	Temperature calibration set

Approval

CE mark	✓
---------	---

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	215×150×30 mm
Net weight	0,2 kg
Shipping method	Parcel service
Net weight approx.	0,20 kg
Gross weight approx.	0,35 kg
Shipping weight	0,35 kg

Pictograms