

KERN MSB 80

KERN

When every centimetre counts – mechanical height rod approved as a medical device for professional use in medical diagnostics



Shipping weight 1,4 kg

Pictograms

STANDARD



1 DAY

Category

Brand	KERN
Product category	Medical scale
Product group	Height rod
Product family	MSB

Measuring System

Units	cm
Division	1 mm
Measuring range length [Min]	10 cm
Measuring range length [Max]	80 cm
Tolerance measuring rod	5 mm

Approval

CE mark	✓
Approved as a medical device	✓

Construction

Dimensions device (W×D×H)	880×330×100 mm
Dimension housing (W×D×H)	882×70×28 mm

Environmental conditions

Use temperature [Min]	10 °C
Use temperature [Max]	40 °C
Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	135×615×85 mm
Net weight	0,6819 kg
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	0,95 kg