

KERN MSC 100

KERN

When every centimetre counts – mechanical height rod approved as a medical device for professional use in medical diagnostics



Shipping weight 0,919 kg

Pictograms

STANDARD



Category

Brand	KERN
Product category	Medical scale
Product group	Height rod
Product family	MSC

Measuring System

Units	inch/ft cm
Division	1 mm
Measuring range length [Min]	7,5 cm
Measuring range length [Max]	100 cm
Tolerance measuring rod	5 mm
Measuring range measuring rod (cm)	7,5-100cm

Approval

CE mark	✓
Approved as a medical device	✓

Construction

Dimensions device (W×D×H)	1040×290×55 mm
Dimension housing (W×D×H)	1020×290×50 mm

Environmental conditions

Use temperature [Min]	10 °C
Use temperature [Max]	40 °C
Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	1085×70×40 mm
Net weight	0,694 kg
Shipping method	Parcel service
Net weight approx.	0,70 kg
Gross weight approx.	0,95 kg