

KERN MSW 200

KERN

Convenient circumference measuring tape with practical pull-out mechanism



Net weight approx.	0,05 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

Pictograms

STANDARD



Category

Brand	KERN
Product category	Medical scale
Product group	Body size measurement
Product family	MSW

Measuring System

Default unit	mm
Units	mm
Division	0,1 cm
Measuring range length [Min]	10 cm
Measuring range length [Max]	205 cm
Tolerance measuring rod	5 mm
Measuring range measuring rod (cm)	0-205cm @ 0,1cm ± 0,5cm

Construction

Dimensions device (W×D×H)	55×90×20 mm
Dimension housing (W×D×H)	55×90×20 mm
Measurement alignment	horizontal

Environmental conditions

Use temperature [Min]	10 °C
Use temperature [Max]	40 °C
Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	70×100×25 mm
Net weight	0,045 kg
Shipping method	Parcel service