

KERN OBB-A 1 156

KERN

5W LED Epi Fluorecence unit (B/G) including centering objective



Category

Brand	Optics
Product category	Optical Instrument
Product group	Fluorescence unit

Power Supply

Supplied power supply	Power supply unit
-----------------------	-------------------

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	395×390×310 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	6 kg
Shipping weight	9,6 kg

Pictograms

STANDARD

