

# KERN OBB-A1157

**KERN**

5W LED Epi Fluorecence unit (B/G) including centering objective



## Category

Brand	Optics
Product category	Optical Instrument
Product group	Fluorescence unit

## Power Supply

Supplied power supply	Power supply unit
-----------------------	-------------------

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	395×390×310 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	6 kg
Shipping weight	9,6 kg

## Pictograms

### STANDARD

