

KERN OBB-A 1466



Blue filter, suitable for models: OBS 104, OBS 106, OBE-1



Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope filter

Construction

Filter type	Standard
Filter diameter	32 mm
Filter colour	blue

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	50×50×2 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,20 kg
Shipping weight	0,2 kg

Pictograms

STANDARD



1 DAY